



Working with the Haemonetics® PCS®2 - Operation Manual -

Printed in France
Haemonetics Corporation
400 Wood Road
Braintree, MA, 02184, USA
©2002, Haemonetics International. All rights reserved.

CE 0123
Authorized European Representative
HAEMONETICS U.K. LTD.
5 Ashley Drive, Bothwell
Scotland G71 8DA
Great Britain

P/N 85266-30, Manual revision: A
April 2002

CONSUMER INFORMATION

Proprietary rights

The contents of this manual are property of the Haemonetics Corporation. Haemonetics® and PCS®2 are registered trademarks of the Haemonetics Corporation. Any information or descriptions contained in this manual may not be reproduced and released to any of the general public, or used in conjunction with any professional instruction without written consent of Haemonetics Corporation, USA. Please direct any written inquiries to the appropriate address:

International Headquarters

Haemonetics S.A.
Signy Center, P.O. Box 262
CH-1274 Signy 2 Switzerland
Tel. +41-22-363-9011
Fax +41-22-363-9054

Corporate Headquarters

Haemonetics Corporation
400 Wood Road
Braintree, MA, 02184, USA
Tel. +1-781-848-7100
Fax +1-781-356-3558

Legal disclaimer

This manual is intended for use as a guide, uniquely for material as supplied by the Haemonetics Corporation. It provides the operator with necessary information to safely carry out specific procedures and satisfactorily maintain Haemonetics-produced equipment. The manual is to be used in conjunction with instruction and training as supplied by qualified Haemonetics personnel.

Haemonetics guarantees its products when correctly used by a properly trained operator. Any failure to respect the procedures as described could result in impaired function of the equipment, as well as in injury to the operator and/or patient /donor. Haemonetics accepts no responsibility for problems resulting from failure to comply with prescriptions as outlined by the company. Any modifications estimated as necessary by the customer should be evaluated by a Haemonetics Clinical Specialist.

Safe utilization of Haemonetics material and equipment requires the operator to correctly handle and dispose of blood-contaminated material. The operator of any Haemonetics equipment must understand and implement the local policies and standard operating procedures concerning the handling of blood-contaminated material, as well as blood products.

It remains solely the responsibility of the customer to fully assess and ensure the safety of any products obtained from Haemonetics prescribed procedures, prior to further application or use. Haemonetics declines any responsibility for choices made by the consumer concerning the utilization of products and by-products.

In addition, it is the responsibility of the apheresis center using Haemonetics equipment and material to inform the donor about the risks involved with any apheresis procedure. Prior to initiating any procedure, the apheresis center is responsible to verify that the donor understands these risks and consents to the procedure.

Haemonetics worldwide locations

Haemonetics Asia Inc.

Taiwan Branch
26F-1, No. 102 Roosevelt Road Sec. 2
Taipei, Taiwan
Tel. +886-2-2369-0722
Fax +886-2-2364-3698

Haemonetics GesmbH

Handelsges.m.b.H.
Berlagasse 45/B2-02
A-1210 Wien, Austria
Tel. +43-1-294-29-00
Fax +43-1-294-29-05

Haemonetics Belgium NV

Leuvensesteenweg 542-BP. 14
Planet II Complex
B-1930 Zaventem, Belgium
Tel. +32-2-720-7484
Fax +32-2-720-7155

Haemonetics BV

C/O CITCO - WTC, PB 7241
Strawinskylaan 1725
1007 JE Amsterdam
The Netherlands
Tel. +31-35-602-3425
Fax +31-35-602-4198

Haemonetics Medical Devices (Shanghai) International Trading Co. Ltd.

Room 28032, Shanghai HSBC Tower
101 Yin Cheng East Road
Shanghai 200120, PRC
Tel. +86-21-506-63366
Fax +86-21-684-13688

Haemonetics CZ, spol. S.r.o

Ptašínského C.8
60200 Brno, Czech Republic
Tel. +42-05-412-122400
Fax +42-05-412-122399

Haemonetics France S.A.R.L.

46 bis, rue Pierre Curie
Z.I. Les Gatines
F-78370 Plaisir, France
Tel. +33-1-30-81-4141
Fax +33-1-30-81-4130

Haemonetics GmbH

Rohrauerstrasse 72
D-81477 München, Germany
Tel. +49-89-785-8070
Fax +49-89-780-9779

Haemonetics Hong Kong Ltd.

Suite 1314, Two Pacific Place
88 Queensway, Hong Kong
Tel. +852-286-89218
Fax +852-280-14380

Haemonetics Italia S.R.L.

Via Donizetti, 30
20020 Lainate (MI), Italy
Tel. +39-2-935-70113
Fax +39-2-935-72132

Haemonetics Japan K.K.

Kyodo Building 3F
16, Ichiban-cho, Chiyoda-ku
Tokyo, Japan, 102-0082
Tel. +81-3-3237-7260
Fax +81-3-3237-7330

Haemonetics Scandinavia AB

Beta Huset, Ideon
Scheelegatan 17
S-223 70 Lund, Sweden
Tel. +46-46-286-2320
Fax +46-46-286-2321

Haemonetics (UK) Ltd.

Beechwood House
Beechwood Estate
Elmete Lane, Roundhay
Leeds LS8 2LQ, United Kingdom
Tel. +44-113-273-7711
Fax +44-113-273-4055

Haemonetics S.A.

Signy Centre
P.O. Box 262
CH-1274 Signy 2, Switzerland
Tel. +41-22-363- 9011
Fax +41-22-363- 9054

Table of Contents

Chapter 1 Explaining General Information

PROVIDING AN OVERVIEW	1-2
What is apheresis technology?	1-2
What is the purpose of this manual?	1-2
What is the Haemonetics Plasma Collection System 2?	1-3
What are the characteristics and features of the PCS2?	1-3
What are the special features of the PCS2?	1-4
What is required to perform a procedure?	1-4
UNDERSTANDING THE USE OF SYMBOLS	1-5
Symbols found in this document	1-5
Symbols found on the device	1-5
Symbols found on disposable packaging	1-8
LISTING DEVICE SPECIFICATIONS	1-9

Chapter 2 Describing the PCS2 Device Components

PRESENTING THE PCS2 DEVICE COMPONENTS	2-3
DESCRIBING THE CENTRIFUGE SYSTEM	2-4
System-sealing mechanism	2-5
Centrifuge well	2-6
Centrifuge base	2-6
DESCRIBING THE PCS2 CABINET COMPONENTS	2-7
Optical line Sensor	2-7
Weigher	2-7
Pumps	2-8
Valves	2-9
Donor flow lights	2-10
Air detectors	2-11
Pressure monitors	2-13
Blood filter holder	2-14
Solution-bag poles (2)	2-14
Power entry module	2-15
Power cord	2-15
Pressure cuff	2-15
Biohazard waste bag	2-15
Communication box/data card (optional)	2-16
Bar code reader (optional)	2-16

DESCRIBING THE PCS2 CONTROL PANEL	2-17
Display screen	2-17
Mode control keys	2-19
Protocol key	2-19
Pump control keys	2-20
Programming keys	2-20
Cuff key	2-21
Valve control keys	2-21

Chapter 3 Maintaining the PCS2 Equipment

CLEANING PROCEDURES	3-2
Cabinet, control panel and valves	3-2
Pressure monitors	3-3
Air detectors	3-3
Optical sensors	3-3
Fluid detector	3-4
Centrifuge components	3-4
Pumps	3-5
Filter screens	3-5
Barcode reader	3-5
CUSTOMER SERVICE	3-6
Clinical training	3-6
Field service	3-6
Returned Goods Authorization system	3-6
HAEMONETICS® CLEANING AND MAINTENANCE RECORD	3-8

Chapter 4 Ensuring Safety and Quality for a PCS2 Procedure

HANDLING THE PCS2 EQUIPMENT	4-2
Storing the PCS2 device and material	4-2
Inspecting the material	4-2
PREVENTING PROBLEMS DURING A PCS2 PROCEDURE	4-3
Understanding the risk of hemolysis	4-3
Avoiding the consequences of flow restriction	4-3
Avoiding bowl misalignment	4-4
Avoiding overheating due to mechanical situations	4-4
Controlling for Red Cell Overrun	4-4
WARNINGS FOR THE OPERATOR	4-5
Electrical shock hazards	4-5
Leakage current control	4-5
Mechanical hazards/rotating parts	4-5
Power outlet connection	4-5
Communicable disease precautions	4-6

Chapter 1

Explaining General Information

PROVIDING AN OVERVIEW	1-2
What is apheresis technology?	1-2
What is the purpose of this manual?	1-2
What is the Haemonetics Plasma Collection System 2?	1-3
What are the characteristics and features of the PCS2?	1-3
What are the special features of the PCS2?	1-4
What is required to perform a procedure?	1-4
UNDERSTANDING THE USE OF SYMBOLS	1-5
Symbols found in this document.	1-5
Symbols found on the device	1-5
Symbols found on disposable packaging	1-8
LISTING DEVICE SPECIFICATIONS	1-9

PROVIDING AN OVERVIEW

What is apheresis technology?

Apheresis is the general term used to describe the separation, selective removal and collection of one or more of the individual components which together form whole blood. This term can be subdivided into two categories:

- **Cytapheresis:** selective removal of one or more of the formed, cellular components of whole blood. These elements include erythrocytes, thrombocytes and leukocytes.
- **Plasmapheresis:** selective removal of plasma, the liquid suspension medium of blood. Plasma contains elements referred to as fractionable components, such as clotting proteins and immunoglobulins.

Apheresis technology permits:

- The collection and separation of whole blood.
- The selective removal and collection of specific components.
- The subsequent return of the non-selected components to the donor or patient.

What is the purpose of this manual?

This manual is intended to supply anyone involved in using Haemonetics equipment with the essential tool for safe and successful operation – *information*. Using this tool of information, the operator can acquire knowledge to be applied throughout all levels of operating experience. This body of information should be consulted whenever necessary, starting from the initial contact with Haemonetics technology to attain:

- An awareness of the purpose of the device and the implications of its collection procedures for the donor and the apheresis center.
- An understanding of how to safely operate the Haemonetics system, correctly install the appropriate disposable material, and troubleshoot any difficulties.
- An ability to consistently apply the principles behind safe operation, proper maintenance and correct handling to ensure optimal, quality apheresis results.

What is the Haemonetics Plasma Collection System 2?

Using updated apheresis technology, Haemonetics has produced the *PCS2* - a compact, lightweight plasmapheresis system which is as easy and safe to use as it is technologically advanced.

The *PCS2* automated apheresis technology provides the operator with a maximum degree of flexibility in any type of plasmapheresis location. The plasma collected may be designated for use in therapeutic transfusion. It may also be conserved, used as source plasma and subsequently fractionated into plasma-derived products.

What are the characteristics and features of the *PCS2*?

The *PCS2* is appropriately called a *Plasma Collection System* because it consists of distinctive *parts* which collectively function as a *whole system* to produce a designated final product:

- The automated plasmapheresis device developed by Haemonetics called the *PCS2*.
- The process designed by Haemonetics to gather plasma from a donor called a *collection procedure*.
- The single-use collection material manufactured by Haemonetics called a *disposable set*.

Once the operator has initiated a *PCS2* procedure, plasma collection will proceed automatically. The appropriate amount of anticoagulant solution will be mixed in the disposable tubing with whole blood from the donor. This anticoagulated blood will be drawn into a disposable collection bowl and separated by centrifugal force into its various components.

When the bowl reaches its collection capacity, the plasma component will exit the bowl and be directed into a plasma collection container for conservation. Non-selected blood components will be returned to the donor. This cycle will be repeated until the desired amount of plasma is collected.

The choice of the disposable collection material will depend on the desired collection product. The *PCS2* technology also provides the operator with the option to infuse saline solution along with the blood components to the donor at different points of a procedure, depending on the type of disposable bowl in use.

Haemonetics has designed the *PCS2* technology with a degree of automation which permits the operator to interact with the device. The operator should remain attentive to the display screen messages while monitoring the status of the donor. It is possible to modify certain aspects of the collection procedures, based on the needs and requirements of the individual donor and the selected material.

What are the special features of the PCS2?

Haemonetics has incorporated advanced technological features into the portable PCS2 design. Examples of these features, which ensure safety for the donor and permit efficient time-management for the operator, are:

- Self-loading pumps.
- Advanced optical sensors.
- Donor-line tubing pressure monitor.
- Communication data box or internal data card.
- Barcode reader.

What is required to perform a procedure?

PCS2 collection procedures are quick and simple to perform. The following material is required to perform a PCS2 procedure:

- A PCS2 disposable set designed for the selected procedure.
- Venipuncture materials and hemostats.
- Appropriate anticoagulant solution.
- 0.9% normal saline (optional).

The operator will need to:

- ➔ Install the appropriate disposable set.
- ➔ Modify any settings if necessary.
- ➔ Perform a single venous puncture, prior to initiating a procedure.

Plasma collection will proceed automatically until the end-collection target has been reached.

UNDERSTANDING THE USE OF SYMBOLS

Symbols found in this document

The terms *note*, *caution* and *warning* are used in this manual with the following symbols to emphasize certain details for the operator.



Note: Provides useful information regarding a procedure or operating technique when using Haemonetics material.



Caution: Advises the operator against initiating an action or creating a situation which could result in damage to equipment, or impair the quality of the by-products; personal injury is unlikely.



Warning: Advises the operator against initiating an action or creating a situation which could result in serious personal injury to either the donor or the operator.

Symbols found on the device

The descriptions of the following symbols are based on information provided in the following documents:

- IEC Standard 60601-1, Medical Electrical Equipment, Part 1: General requirements for safety.
- IEC Standard 60417-1, Graphical symbols for use on equipment, Part 1: Overview and application.



Type BF applied part

This symbol indicates that the applied portion (i.e. the part which comes in contact with the donor) of the device is electrically isolated. The device has an internal electrical power source providing adequate protection against electrical shock, in particular pertaining to acceptable leakage current and the reliability of the protective earth connection.




Protective earth (ground)

Used to identify any terminal intended for connection to an external conductor, for protection against electrical shock in case of a fault.

**Alternating current**

Used to indicate on the rating plate that the device is suitable for alternating current only.

**Fuse symbol**

Used to identify fuse boxes or the location of a fuse box.

**Power OFF**

Position of the main power switch indicating disconnection from the mains.

**Power ON**

Position of the main power switch indicating connection to the mains.

IPX1**Protection against ingress of liquid**

Indicates that the enclosure of the device is designed to provide a specified degree of protection against harmful ingress of water or liquid into the equipment (under applicable conditions).

**Attention** (Consult accompanying documents)**Non-ionizing electromagnetic radiation**

Used to specify RF transmission for data communication.

The following symbols have been designed for devices manufactured by Haemonetics:



Bar-code reader connection



RS232 connection



RS232 connection with power to one pin



Pressure cuff connection

Symbols found on disposable packaging


The following symbols are used by Haemonetics on disposable set packaging.

REF

CATALOG NUMBER



EXPIRATION DATE



Lot Number



Sterilized by exposure to Ethylene Oxide



Fluid path STERILE by exposure to Ethylene Oxide



Sterilized by exposure to Gamma irradiation



Fluid path STERILE by exposure to Gamma irradiation




DO NOT REUSE



Caution: consult operator manual for instructions



Storage conditions, humidity level



Storage conditions, temperature level

LISTING DEVICE SPECIFICATIONS

The approximate weight and dimensions of the PCS2 device are as follows:

Characteristics	Values	
	Cabinet cover open	Cabinet cover closed
Height	63 cm	44 cm
Width	55 cm	
Depth	55 cm	32 cm
Depth with communication box	55 cm	37 cm
Weight	26.4 kg	
Weight with communication box	27.4 kg	

The following environmental conditions should be respected pertaining to operation and storage of the PCS2 device:

Conditions	Values
Ambient operating temperature	+18° C to +27° C
Tested storage temperature	0° C to + 40° C
Storage humidity level	Maximum relative humidity rate of 90%, <i>non-condensing</i>

The electrical specifications for operating the PCS2 device are as follows:

Characteristics	Values (relative to input voltage)	
Input voltage	230 VAC \pm 10%	110 VAC \pm 10%
Operating current	~1.9 A	~ 2.6 A
Fuse rating	F2.5 A @ 250 V	F5.0 A @ 250 V
Operating frequency range	50 - 60 Hz	50 - 60 Hz
Maximum leakage current	500 μ A	100 μ A



Note: Haemonetics will regulate the proper voltage setting upon installation. The power source used must be properly grounded.



Caution: *The PCS2 device must be operated in an environment compatible to the requirements of the IEC 60601-1-2 Standard, Electromagnetic compatibility.*

Mobile RF communication equipment not approved by Haemonetics and portable communication equipment can affect the PCS2 device. Any accessories and cables not approved by Haemonetics used in conjunction with the device may increase hazards and influence compatibility with EMC requirements. Therefore, non-approved accessories and cables must not be used.

In addition, the PCS2 device and accessories must not be placed directly adjacent to, or top of other equipment, unless specifically approved by Haemonetics.

Chapter 2

Describing the PCS2 Device Components

PRESENTING THE PCS2 DEVICE COMPONENTS.....	2-3
DESCRIBING THE CENTRIFUGE SYSTEM.....	2-4
System-sealing mechanism.....	2-5
Centrifuge well.....	2-6
Centrifuge base.....	2-6
DESCRIBING THE PCS2 CABINET COMPONENTS.....	2-7
Optical line Sensor.....	2-7
Weigher.....	2-7
Pumps.....	2-8
Valves.....	2-9
Donor flow lights.....	2-10
Air detectors.....	2-11
Pressure monitors.....	2-13
Blood filter holder.....	2-14
Solution-bag poles (2).....	2-14
Power entry module.....	2-15
Power cord.....	2-15
Pressure cuff.....	2-15
Biohazard waste bag.....	2-15
Communication box/data card (optional).....	2-16
Bar code reader (optional).....	2-16
DESCRIBING THE PCS2 CONTROL PANEL.....	2-17
Display screen.....	2-17
Mode control keys.....	2-19
Protocol key.....	2-19
Pump control keys.....	2-20
Programming keys.....	2-20
Cuff key.....	2-21
Valve control keys.....	2-21

A. Cabinet

1. Centrifuge
2. Line sensor
3. Weigher
4. Anticoagulant (AC) pump
5. Blood pump
6. Donor valve (red)
7. Plasma valve (yellow)
8. Saline valve (white)
9. Donor flow lights (x2)
10. Anticoagulant line air detector (ACAD)
11. Blood line air detector (BLAD)
12. Donor line air detector 1 (DLAD1)
13. Donor line air detector 2 (DLAD2)
14. Blood filter holder
15. AC solution pole
16. Saline solution pole
17. Donor pressure monitor (DPM)
18. System pressure monitor (SPM)

B. Control Panel

19. Display screen
20. Mode control keys
21. Protocol key
22. Pump control keys
23. Programming keys
24. Cuff key
25. Valve control keys




Figure 2-1, PCS2 device components

PRESENTING THE PCS2 DEVICE COMPONENTS

The components of the PCS2 device will be presented in this chapter according to where they are located on the device:

- The centrifuge system.
- The cabinet components.
- The control panel.

1. Centrifuge system
2. Line sensor
3. Weigher
4. AC pump
5. AC pump tubing guide
6. Blood pump
7. Blood pump tubing guide
8. Donor valve (red)
9. Plasma valve (yellow)
10. Saline valve (white)
11. Donor flow lights (x2)
12. ACAD (Anticoagulant air detector).
13. BLAD (Blood line air detector).
14. DPM (Donor pressure monitor).
15. SPM (System pressure monitor).
16. DLAD X 2 (Donor line air detectors).
17. Blood filter brackets
18. Tubing guide
19. AC solution pole
20. Saline solution pole




Figure 2-2, PCS2 top deck view

The disposable set elements will be explained in greater detail in the chapter “Describing the PCS2 Disposable Collection Material”.

As an explanation for the references made to the disposable tubing in this chapter:

Donor-line tubing refers to the tubing which is either:

- Transporting blood away from the donor before entering the bowl.
- Transporting blood from the bowl before return to the donor.

Effluent tubing refers to the tubing exiting the bowl in the direction of the plasma collection container.

DESCRIBING THE CENTRIFUGE SYSTEM

The centrifuge system of the PCS2 device is designed to hold a disposable bowl in which the blood components can be spun from a range of 3000 to 8000 revolutions per minute. This centrifugal force will separate anticoagulated whole blood into its various components.

The PCS2 centrifuge system consists of:

- A system-sealing mechanism.
- The centrifuge well.
- The centrifuge base.

1. Split hinged lid
2. Fluid detector
3. Optical bowl sensor
4. Locking knob
5. Centrifuge well
6. Centrifuge chuck
7. Centrifuge base




Figure 2-3, PCS2 centrifuge system components

System-sealing mechanism

The PCS2 centrifuge contains a split, hinged lid (or *cover*) and a locking knob. These components “seal” the system by:

- Securing the contact of the disposable bowl with the centrifuge base.
- Isolating the spinning bowl from the operator.

Centrifuge cover

The centrifuge lid, referred to as the *cover*, has tabs located on the rimmed portion of each split side. The split halves are attached to the centrifuge rim by a hinge. As the halves of the lid are lowered to meet the rim, the tabs must be firmly pressed together in order to completely close the lid and provide a seal around the stationary head of the disposable bowl.

The split halves of the lid are made from a durable, transparent material, allowing the operator to observe changes in the bowl contents as the centrifuge spins.

Locking knob

The knob is positioned on the rim of the centrifuge well.

Once the lid has been fully closed, the knob requires a series of turns in a clockwise direction to lock the centrifuge and thus completely seal the system.

To unlock the centrifuge, the operator should turn the knob in a counter-clockwise direction until the split halves can be separated, then lifted, to open the lid.



Figure 2-4, PCS2 centrifuge cover

Centrifuge well

The PCS2 centrifuge well is designed with the following components:

Optical bowl sensor

There is an optical sensor located on the upper portion of the centrifuge well. The sensor is aimed at the core of the bowl and will measure optical reflection as the various blood components pass in front of the optical beam.



Note: The interface between the optical sensor in the centrifuge well and the contents of the bowl is often referred to as the “bowl optics readings” and will be discussed further in the chapter “Understanding a PCS2 Collection Procedure”.

Fluid detector

The PCS2 centrifuge well is equipped with an electronic fluid detection system designed to detect the presence of liquid. The detector is mounted on the wall of the centrifuge well. The PCS2 safety system will automatically stop the centrifuge and the pumps if there is contact between liquid of any sort and the fluid detector.

1. Optical bowl sensor
2. Fluid detector
3. Centrifuge chuck




Figure 2-5, PCS2 centrifuge well components

Centrifuge base

The centrifuge base contains a centrifuge chuck designed to hold the disposable bowl in place during operation.

When installing a bowl, the operator should exert a downward pressure on the head of the bowl and ensure that the bowl is completely seated. A suction force will be created, between the base of the bowl and the chuck, to hold the bowl in place. The bowl will be completely secured once the operator has locked the centrifuge lid.

To remove a bowl at the end of a procedure, the operator should open the centrifuge lid and pull upward on the bowl.



Warning: The PCS2 device is equipped with a safety feature which will not allow the centrifuge to spin if the lid has been improperly closed. It is unlikely that a properly installed centrifuge bowl will become unaligned as it spins.

If the operator does notice anything unusual, under no circumstances, should the operator attempt to open the centrifuge lid if the bowl is still spinning. The operator must ensure that the bowl has come to a complete stop before attempting to open the centrifuge for any reason.

DESCRIBING THE PCS2 CABINET COMPONENTS

Optical line Sensor

Located on the right side of the PCS2 top deck is the optical line sensor which monitors the blood components passing through the effluent tubing.



Caution *The line sensor will not provide accurate readings if the optical lens it is obstructed in any way; thus the lens must be cleared of any extraneous substances to ensure proper functioning of the system.*

Weigher

The *weigher* is the term used by Haemonetics to describe the PCS2 component which measures in grams the contents of the plasma collection container(s) placed on the *weigher arm*. When the Draw key is pressed to initiate a procedure, the weigher will automatically tare, or deduct the weight of the empty plasma collection container. Thus, the weight of the container will not be included in the weight displayed for the PCS2 procedure statistics.

To ensure optimal accuracy from the weigher during a collection procedure:

- The weigher arm must be fully extended, positioned at a 90 degree angle to the PCS2 top deck, prior to the system self-tests.
- The plasma collection container must hang freely without any contact with the PCS2 cabinet.



Caution: *The operator must be careful to not touch the weigher once the weight of the plasma container has been tared. This could affect the programmed target plasma quantity and could result in an excessive collection of plasma during a Draw cycle. If the weigher senses a decrease in the weight during a Return cycle, an error message will be displayed.*

To ensure that the weigher arm is within the appropriate calibration range, the operator has the option to verify the weight displayed during the READY mode, prior to the first Draw cycle, as follows:

- ➔ Remove the plasma container from the weigher arm.
- ➔ Ensure that the weigher arm is fully extended at a 90 degree angle to the top deck.
- ➔ Hang a certified weight (*not exceeding 1300 grams*) from the empty weigher arm.
- ➔ Note the value displayed (*and convert to grams if displayed in ml*).
- ➔ Replace the plasma container on the correctly positioned weigher arm.

The measurement of the weight displayed in grams should be within 1% of the certified weight.

Pumps

Located on the left side of the PCS2 top deck are two pumps which use peristaltic movements to displace fluids through the disposable-set tubing.




Figure 2-6, PCS2 pump rotor

Anticoagulant pump

The AC pump, designated by the color blue, moves AC solution between the AC solution bag and the needle connector of the donor line tubing.

Blood pump

The Blood pump, designated by the color red, moves fluids between the donor and the centrifuge bowl.

The pumps function during the different modes of PCS2 operation as follows:

- When loading the disposable tubing, the AC pump and the Blood pump turn simultaneously to thread the tubing onto the pump rotors.
- During the PRIME mode, the AC pump and the Blood pump turn simultaneously to provide the inlet side of the donor line tubing with AC solution.
- During the DRAW mode, the AC pump and the Blood pump turn simultaneously. The Blood pump pulls anticoagulated whole blood through the blood filter of the disposable set into the centrifuge bowl.
- During the RETURN mode, the Blood pump pulls the remaining blood components from the centrifuge bowl and re-infuses them to the donor.



Note: The two pumps will rotate at different speeds during the PCS2 modes of operation, depending on the AC/Blood pump ratio procedure parameter setting.

Pump tubing guide

Next to each of the PCS2 pumps is a tubing guide. It will secure the disposable tubing during pump autoloading, as well as during the collection procedure.

Valves

There are three valves located on the PCS2 top deck which automatically control the flow of fluids through the disposable set tubing. The valves are color-coded according to their specific functions.

The PCS2 safety system will control the valves during the self-diagnostic tests. Once the operator has selected a collection procedure, the appropriate valves will automatically open, in preparation for loading the disposable set tubing. If the disposable tubing needs an adjustment during a procedure, it is possible to open a valve manually by pressing the valve lever located at the top of each valve.



Warning: Any manual adjustment to a valve should be attempted only if the PCS2 device is POWERED-OFF, in the READY mode, or when the pumps are stopped. At any other time, the PCS2 safety system will be alerted and will interrupt the procedure. Manipulating a valve could lead to flow problems, and eventually cause hemolysis.

Donor valve (red)

The donor valve is located the furthest to the left on the PCS2 top deck.

- During the DRAW mode it remains open so that anticoagulated whole blood can flow from the donor into the centrifuge bowl.
- During the RETURN mode it remains open so that components remaining in the disposable tubing and centrifuge bowl can be returned to the donor.

Plasma valve (yellow)

The plasma valve is located between the donor valve and the saline valve.

- During the DRAW mode it remains open to direct plasma and air flowing through the effluent tubing into the plasma collection container.
- During the RETURN mode it remains open, except for a brief period when the effluent line is cleared during a standard bowl procedure. It is closed during automatic saline compensation to the donor.

Saline valve (white)

The saline valve is located the furthest to the right on the PCS2 top deck. It will be used to direct saline solution through the tubing, with the exception of a multi-bag plasma collection procedure.

- During the DRAW mode it will be closed.
- During the RETURN mode it will be open to permit the passage of fluid through the effluent tubing.



Note: When the selected PCS2 procedure uses more than one plasma collection bags, plasma will be distributed equally between the bags, using the plasma valve to route half of the plasma being collected, and the saline valve to route the other half. The two valves will not be open at the same time during DRAW, therefore automatic saline compensation is not a possibility when more than one plasma bag is used.

Donor flow lights

These color-coded lights, located on both sides of the PCS2 top deck, indicate donor blood-flow status during the DRAW and RETURN modes. They are contained in a rectangular panel on the PCS2 top deck.




Figure 2-7, PCS2 donor flow lights

DRAW mode

- The GREEN light indicates that donor blood flow is adequate for the blood pump to maintain an adequate speed.
- The YELLOW light indicates that donor flow is decreasing to a rate of less than 2/3 of what is required to maintain an adequate speed.
- The RED light indicates that blood is not flowing adequately or not flowing at all from the donor.

If the RED flow light is lit, the Blood pump will automatically stop. The centrifuge will continue to spin to ensure continued separation of the collected blood components. When donor blood flow is restored, the pump will automatically restart.

RETURN mode

- The YELLOW light indicates that the non-selected blood components are being returned to the donor. No other donor light will be visible at this time.



Note: If any of the DRAW mode lights are lit, the donor can promote blood flow by clenching and relaxing the hand below the needle site. When the RETURN yellow light is lit, the donor should NOT do this, because the blood components in the bowl are being returned. The operator should instruct the donor to observe the differences in the lights and act accordingly.

Air detectors

The PCS2 is equipped with an assembly of ultrasonic sensors designed to detect the presence of air, bubbles or foam in the fluids flowing through the disposable set tubing. If air is detected outside of the normal range during any mode (PRIME, DRAW OR RETURN), the detectors will:

- Alert the PCS2 safety system.
- Stop the operation in progress.
- Provide the operator with an error message and an audible alarm.




Figure 2-8, Example of a PCS2 air detector

Anticoagulant Air Detector (ACAD)

The passage of anticoagulant solution from the bag into the PCS2 system occurs over a series of steps. The ACAD will monitor the AC solution line throughout the entire procedure.

The ACAD is located on the top deck of the PCS2 cabinet adjacent to the AC pump and will function once the AC solution has passed through the AC pump tubing.

Blood Line Air Detector (BLAD)

This air detector, located on the top deck of the PCS2 cabinet to the right of the blood pump, will serve a dual purpose to the operator. The BLAD will remain active throughout the entire procedure, but will provide a specific function at the following moments:

- During a Draw cycle, the BLAD will detect the presence of fluid passing through the blood line tubing. This allows the system to account for the volume of blood being pumped.
- During a Return cycle, the BLAD will note the presence of any air in the tubing leaving the centrifuge bowl. This line contains the blood being returned to the donor and will pass through the donor valve after the BLAD. When the BLAD has detected air in the tubing within normal limits, this will signal that the bowl is empty and the Return cycle will be terminated.



Warning: Air detected (or lack of air detection) by the BLAD, outside of normal limits, will stop the collection procedure and alert the operator.

Donor line air detectors (DLAD1 and DLAD2)

The two donor line air detectors are located on the left side of the PCS2 front panel; the DLAD1 is located above the DLAD2. Both air detectors monitor the donor line between the donor and the disposable set blood filter as follows:

- During the PRIME mode, the pumps draw AC solution into the donor line up to the DLAD. When the DLAD note fluid (or an absence of air), the detectors signal to PCS2 software that the line has been primed with anticoagulant solution and is prepared for the Draw cycle.
- During the DRAW mode, the DLAD will monitor the tubing containing anticoagulated whole blood and alert the operator if any air has been introduced into the system.
- During the RETURN mode, the DLAD will monitor the donor line as it carries blood components being returned to the donor. The DLAD monitor the line for any air which may have passed into the system undetected (*probability very low*) by the BLAD.



Warning: In the case of any air detection alarm, the operator must respond immediately, note the source and take appropriate action. Further information is provided in the chapter “Troubleshooting for a PCS2 Procedure.”

During a Return cycle, if either the DLAD1 and/or DLAD2 produce an air detection alarm, this could indicate a failure of the BLAD. The operator should carefully note the source of air detected; no blood should be sent to the donor until all air bubbles have been removed from the line.

Haemonetics recommends the following operator actions to remove any air bubbles detected in the tubing between the BLAD, DLAD1 and the DLAD2:

- ➔ Press the Draw key until blood enters the bowl, to send any air bubbles to the bowl.
- ➔ Continue with a Return cycle only after any air bubbles have been removed.

1. DLAD1
2. DLAD2
3. Blood filter brackets
4. AC solution pole
5. DPM
6. SPM




Figure 2-9, PCS2 front/side view

Pressure monitors

The electronically controlled pressure monitors function with the correlating filter on the disposable set to measure pressure in the disposable tubing. The pressure monitors provide feedback to the system about the flow of blood components. The PCS2 programming will automatically regulate the speed of the pumps based on this information.

Donor pressure monitor (DPM)

The DPM, located on the left side of the PCS2 top deck, measures pressure in the donor line tubing. The information is depicted on the display screen using a bar graph. The bar graph is visible on the screen when donor pressure is adequate to maintain the programmed pump speed. The bar graph will not be visible if donor-line pressure is below what is required to maintain the programmed pump speed.

Variations will exist in the readings, depending on the operating mode. The PCS2 software is programmed to detect a range of normal values. If a pressure reading varies outside of this range, the PCS2 safety system will stop the pumps and provide an explanatory screen message, as well as an intermittent alarm.



Caution: *Once the DPM and the disposable set filter have been connected, the system is ready to measure the pressure in the donor line. This connection should not be disrupted at any point. In the event of a power failure, the operator should refer to the chapter “Troubleshooting During a PCS2 Procedure” for details about handling the DPM and the donor line tubing. The DPM tubing must be clamped before removing the filter to ensure that the procedure can be recovered.*

DRAW mode

The pressure readings will vary as blood is drawn from the donor. If a significant pressure decrease is detected and the DPM readings drop below a programmed value, the pump speed will automatically decrease until a sufficient pressure increase is measured.

If the donor-line pressure is measured as insufficient, the pumps will stop, the NO FLOW indicator lights will be visible and an explanatory screen message will appear with an alarm. Once pressure is measured to be within normal operating range, the pumps will resume the programmed speed.

RETURN mode

The pressure readings will vary as blood is returned to the donor. If a significant pressure increase is detected, and the pressure readings rise above a programmed value, the pump speed will automatically decrease until a sufficient pressure change is measured. If pressure readings remain high, the Blood pump will stop, and an explanatory screen message will appear with an alarm.

When pressure is measured to be within the normal operating range, the pumps will resume operation until reaching the programmed pump speed.



Warning: The operator must remain aware of the fact that a high pressure warning can indicate a possible flow obstruction and could cause red blood cell hemolysis, and/or damage the vein.

Corrective action is necessary and the operator should consult the chapter “Troubleshooting During a PCS2 Procedure”, as well as the chapter “Ensuring Safety and Quality for a PCS2 Procedure” for further information.

System pressure monitor (SPM) (optional)

The SPM, located on the right side of the PCS2 top deck, measures pressure in the effluent tubing. This measurement verifies that the sterile seal, between the head and the body of the centrifuge bowl, remains intact.

If the SPM detects that pressure in the system increases or decreases abnormally, the PCS2 safety system will stop the pumps and provide an explanatory message with an alarm. Centrifuge function will remain unaffected.

Blood filter holder

Located on the left side of the PCS2 front panel are two brackets designed to secure the blood filter chamber of the disposable set during the procedure.

Solution-bag poles (2)

Located on either side of the PCS2 cabinet is a height-adjustable pole. These poles are used to hang the solution bags during the procedure. The left pole should be used to hang the anticoagulant solution, whereas the right pole should be used to hang the saline solution (if selected).

There is a knob located on the base of each pole. Pulling the knob outward will disengage the contact pin from the pole and allow the operator to extend or retract the pole.

Power entry module

The power entry module is located on the left panel of the device. Externally, the module consists of an ON/OFF power switch and a power-input receptacle for the power cord. Internally, the module contains the fuse panel. It will interrupt power supply to the system in the event of an electrical current surcharge.

The design of the power entry module also provides a filter-effect for the PCS2 device against the effects of a power surge

In the case of an emergency, the ON/OFF switch can be used to stop all device function.




Figure 2-10, Power entry module (PEM)

Power cord

The power cord provided is designed to connect the PCS2 device with an external power source via the power-input receptacle, located on the power-entry module on the left side panel.

Pressure cuff

The tourniquet-style pressure cuff is used to maintain an optimal venous blood-flow from the donor during specific phases of the collection procedure. The cuff should be attached to the PCS2 cuff connector located on the PCS2 rear panel.

Biohazard waste bag

The biohazard waste bag is designed to collect any biologically contaminated material from the centrifuge well in the rare case of a spill or leak. Two biohazard waste bags are supplied with the delivery of each PCS2 device.

A bag must be attached at all times to the centrifuge drain tube, located at the rear of the device. The bag must hang freely, with the clamp open, visible to the operator.



Warning: The biohazard waste bags are not to be used to collect or store apheresis products. When a bag contains evacuated waste products, it must be clamped, removed and properly disposed of, according to the local policies concerning biologically contaminated material. A new bag must be placed before resuming operation.

1. Power entry module
2. Power cord
3. Pressure cuff
4. Biohazard waste bag



Figure 2-11, PCS2 rear panel view

Communication box/data card (optional)

The external communication box or internal data card transfer data from the PCS2 device to another external device such as a printer, or to *HaemoNet*, the Haemonetics communication network. HaemoNet provides any establishment using Haemonetics equipment with the possibility of linking several Haemonetics apheresis devices to a central monitoring computer. Using HaemoNet, procedure data can be exchanged and stored in a database and/or viewed directly.



Note: The PCS2 communication box has been tested according to standards required by EN 60601-1-2 (EMC of medical electrical equipment). The measured error rate of data communicated to HaemoNet at certain specific electromagnetic frequencies rises above the standards. However, there is no impact on the integrity of the procedure information stored in the database. The HaemoNet communication is designed with CRC error checking, performed upon the reception of all data.

Bar code reader (optional)

The barcode reader (when installed), is located on the left panel of the PCS2 and can be used to enter the following types of data directly into the PCS2 data storage memory.

- Disposable set lot and list numbers.
- Anticoagulant and saline bag solution codes.
- Donation number, donor number and operator identification code.

DESCRIBING THE PCS2 CONTROL PANEL

The control panel, located on the inside of the hinged PCS2 cabinet cover, consists of a **display screen** and a **keypad** comprised of several groups of keys. The control panel allows the operator to interact with the system by entering appropriate data and observing feedback. There is a protective plastic coating on the keypad, which allows for efficient cleaning and disinfecting.

1. Display screen
2. Mode control keys
3. Protocol key
4. Pump control keys
5. Programming keys
6. Cuff key
7. Valve control keys




Figure 2-12, The PCS2 control panel components

Display screen

The control panel display screen will provide data related to procedure status, as well as donor status, throughout the entire procedure. The information will be displayed in specific areas as follows:

1. Active state message.
2. Procedure acronym.
3. Operating state icon.
4. Procedure statistics area.
5. DPM bar graph reading.




Figure 2-13, Example of a PCS2 display screen message

Message area

The center portion of the screen display will provide textual information about the current state (or mode) of operation, while the actual cycle in progress will be (visible in the upper right corner).







Procedure acronym

The upper left portion of the display screen will contain the acronym selected to describe the type of procedure in progress.

Display screen icons

These symbols, located on the upper left side of the display screen, provide a pictorial representation of the operating state, or mode, in progress.

Table 2-1: PCS2 display screen icons

Display screen icon	Explanation	Operating state/mode
	Displayed during AC solution priming sequence.	PRIME
	Displayed when the device in a non-active state, or ready for an operator command.	READY
	Displayed as donor blood is being drawn into the centrifuge bowl.	DRAW
	Displayed as blood components/fluids are being returned to the donor.	RETURN
	Displayed briefly at the end of the procedure, prior to the completed procedure statistics displayed.	PROCEDURE COMPLETE
	Displayed when the centrifuge is stopping.	

DPM bar graph

This is a visual representation of the pressure reading in the donor-line tubing, as it measured by the donor pressure monitor. The DPM bar graph (when visible) is located on the right side of the screen. The contrast between filled and non-filled area in the column will vary to depict the fluctuations in the DPM readings.

Procedure statistics area

The lower portion of the screen communicates data to the operator concerning specific measurements and calculations made by the system during a PCS2 collection procedure. These statistics are updated throughout a procedure and concern the following component functions:

- Blood pump speed during the DRAW and RETURN modes.
- Amount of the plasma contained in the collection container.
- Cuff pressure registered for the donor (mmHg).
- Saline solution volume infused (if selected) during the procedure.
- Volume processed (ml).

Mode control keys

These keys, located directly below the PCS2 display screen, are used to regulate the operating state (or mode) of a PCS2 collection procedure.

PRIME key

This key is used to initiate the PRIME mode. The PRIME mode will bring anticoagulant solution from the anticoagulant line tubing into the donor- line tubing.

DRAW key

This key is used to initiate the DRAW mode. The DRAW mode will move anti-coagulated whole blood from the donor through the donor-line tubing into the centrifuge bowl, where plasmapheresis will be initiated.

RETURN key

This key is used to initiate, or resume, the RETURN mode. During automated procedure functioning, a Return cycle is automatically initiated. However this key can be used to provide an early return of the bowl contents to the donor, if necessary.

STOP key

This key is used to immediately stop the centrifuge and pumps.



Caution: *If the STOP KEY has been used during a Draw cycle, the bowl contents should be **returned** to the donor before resuming the collection procedure. Stopping the procedure during the DRAW mode could affect the separation of the blood components in the bowl. This could eventually interfere with the quality of the final collection product (and/or the collection procedure).*

Protocol key

This key is used to select certain PCS2 procedures and any options. Further information is provided in the chapter "Performing a PCS2 Collection Procedure".

Pump control keys

These keys can be used by the operator to manually change the preset pump speed during a collection procedure.

Pump arrow keys

These keys can be used to temporarily modify the default parameter settings by respectively increasing (arrow up) or decreasing (arrow down) the speed in which the pumps will rotate. The adjustment should be made based on individual donor needs during a specific collection procedure.



Caution: *During a Draw cycle, the operator should observe the DPM bar graph and flow indicator lights in order to correctly assess for low donor blood flow. However, before using the arrow keys to adjust for low donor flow, it is important that the operator allow the PCS2 device to first reach the preset target pump speed.*



Warning: **During a Return cycle, if the pump speed is manually decreased, the operator must carefully monitor the venipuncture site, to avoid possible consequences of an infiltrated vein for the donor, such as hematoma.**

Pump start/stop key

This key can be used to either stop the pumps, or re-start the pumps if stopped by the operator.



Caution: *If the pumps have been stopped using this key, and remain stopped for longer than two to three minutes during a Draw cycle, the bowl may become over-packed with red cells because the centrifuge will continue to spin. This can create a potential flow problem during the Return cycle. In this case, the operator should return the bowl contents to the donor before proceeding with a Draw cycle.*

Programming keys

This section of the PCS2 keypad consists of four keys which enable the operator to modify specific PCS2 procedure parameters. Certain system operating parameters have been selected by Haemonetics as default values. These parameters provide optimal results in PCS2 plasmapheresis procedures with the average donor, as well as for average collection requirements.

However, it is possible to alter and subsequently retain the altered parameters for specific collection requirements. Once the program parameters have been consulted and/or modified, the operator can return to the screen depicting the current mode of operation without interruption to the collection procedure. Further information is provided in the chapter "Understanding a PCS2 Collection Procedure."

Modify Program key



This key is used to access and scroll the list of procedure parameters; it can be pressed during any of the operating modes. Each time that this key is pressed, a different program parameter will be displayed on the screen, along with the current setting for that parameter. To modify the parameter setting, the operator should use the Yes/No arrow keys.

YES/NO Arrow keys

These keys have a dual function:

- Provide a response to a question-prompt from the PCS2 software.
- Modify the parameter setting displayed on the screen.

To modify a displayed parameter value, the operator should:

- ➔ Press the Yes  key (arrow up) to increase the value.
- ➔ Press the No  key (arrow down) to decrease the value.

Save Program key

The operator can use this key to retain each modified value in the PCS2 memory. If this key is pressed after a modification, the selected value will become the new system default value until any further modification is made during subsequent PCS2 collection procedures.

Cuff key

This key is located on the keypad next to the saline and plasma valve keys. During procedure operation, the cuff will automatically inflate during the DRAW mode and deflate during the RETURN mode. The cuff cannot be inflated during the RETURN mode. The operator can use the key to manually control the pressure cuff:

- Prior to a procedure when performing the venipuncture
- During the READY or DRAW modes, to modify cuff pressure.

Valve control keys

These keys can be used to manually control what would normally be the automatic control of a valve, as in the case of manual saline compensation before the final Return cycle. When a valve control key has been pressed, the key will be lit, indicating that the valve is open. When it is re-pressed, the valve will close.

Saline key

This key can be used to manually control the saline (white) valve.

Plasma key

This key can be used to manually control the plasma (yellow) valve.

Chapter 3

Maintaining the PCS2 Equipment

CLEANING PROCEDURES	3-2
Cabinet, control panel and valves.	3-2
Pressure monitors	3-3
Air detectors	3-3
Optical sensors	3-3
Fluid detector	3-4
Centrifuge components.	3-4
Pumps	3-5
Filter screens	3-5
Barcode reader.	3-5
CUSTOMER SERVICE	3-6
Clinical training	3-6
Field service	3-6
Returned Goods Authorization system	3-6
HAEMONETICS® CLEANING AND MAINTENANCE RECORD	3-8

CLEANING PROCEDURES

The PCS2 device has been designed to require minimal maintenance for the operator. To maintain the precision function of the PCS2 device, the operator needs to primarily perform routine cleaning procedures of certain key components. A record of routine cleaning should be kept along with any routine or preventive service maintenance performed by a Haemonetics representative, and a form is provided at the end of this chapter.

The frequency of cleaning each individual PCS2 device will depend on the number of procedures performed. Special cleaning needs may arise and should be dealt with promptly. Haemonetics recommends the following routine cleaning schedule for each PCS2 device, **based on an average of three collection procedures per day, or approximately sixty per month.**

- **Daily:** Clean the exterior surfaces as well as the pressure monitors.
- **Weekly:** Clean the air detectors, the optical sensors (line sensor and optical bowl sensor), the fluid detector and the inside of the centrifuge well.
- **Monthly:** Clean the pump rotors and the pump wells.
- **Quarterly:** Clean the filter screens.



Warning: To eliminate the potential danger of electrical shock, the operator must clean the PCS2 device only when it is disconnected from an external power source.

The following list describes the basic material required for routine cleaning.

- Disinfectant cleaning solution, specific for blood born pathogens and compatible for cleaning Lexan® plastic.
- Warm water.
- 70% Isopropyl alcohol.
- Lint-free gauze or cloth (for cleaning and drying).
- Cotton swabs.
- Protective gloves.
- Hexagonal-head wrench #10 (provided with the device).
- Silicon lubricant (used for the O-ring gasket of the centrifuge chuck).
- Phillips-head screwdriver.

Cabinet, control panel and valves

The exterior cabinet, keypad, display screen and valves should be wiped daily, as well as following any spill, using an appropriate cleaning solution.



Caution: Certain cleaning solutions can degrade the Lexan plastic parts of the PCS2 valves and only compatible cleaning solutions should be used.

Pressure monitors

The pressure monitors (DPM/SPM) should be cleaned daily in the following manner:

- Depress and hold the white ring as if installing the disposable filter.
- Wipe the silver rod thoroughly, using a circular motion and warm water.
- Dry the rod and release the pressure on the ring.



Caution: *It is very important to use **only** water on the pressure monitor rod. Alcohol or cleaning solution residue could cause a reaction with the plastic material of the corresponding disposable set filter and affect the function of the filter.*

Air detectors

The DLAD1, DLAD2, and BLAD are designed with a groove to hold the disposable tubing. The contents of the tubing are monitored by the sensors which are located internally on either side of this groove. If a procedure is interrupted due to an air detector alarm, the operator should remove the tubing and clean the groove before continuing the collection procedure.

The operator should use **only** warm water and lint-free gauze to clean and dry inside of the tubing groove. The groove should be kept free of any particles, such as powder residue from disposable gloves, since this could lead to an erroneous detection of air.



Caution: *Alcohol may cause the plastic housing to degrade and must not be used to clean the air detectors.*

Optical sensors

The lenses of the optical sensors must be kept completely free of particles or debris which could produce inaccurate readings and influence the PCS2 device performance. The operator should use **only** water and lint-free gauze to clean and dry the lenses.



Caution: *If cleaning solution should come into contact with the optical sensor lenses, the lenses should be cleaned immediately with lint-free gauze and warm water, then thoroughly dried. Dried cleaning solution can leave an opaque film on the lens.*

Line sensor

The line sensor, located on the PCS2 top deck, contains two very small lenses which are centered on either side of the disposable tubing groove. The operator should carefully pass the gauze through this groove to clean and dry the lenses.

Optical bowl sensor

The optical bowl sensor lens is located in the upper portion of the centrifuge well. The operator should ensure that no spots remain on the lens cover after it has been cleaned and dried.

Fluid detector

The fluid detector is located inside of the centrifuge well. The surface of the detector should be cleaned using a cotton swab moistened with 70% alcohol.

Centrifuge components

Except for the optical sensor and fluid detector, the other centrifuge components can be wiped routinely using the cleaning solution and a lint-free cloth. This includes the centrifuge well, chuck, hinged lid and locking knob.

Haemonetics Technical Services provides silicon lubricant for the O-ring gasket, located at the base of the centrifuge chuck. After a major cleaning, the operator should apply a small amount of the lubricant to the gasket to prevent it from cracking. It is not necessary to remove the gasket when applying the lubricant.

If the case of a fluid spill, the operator should:

- ➔ Power off the device and disconnect it from the external power source before cleaning.
- ➔ Ensure that the biohazard waste bag is attached to the drain tube and hangs freely.
- ➔ Wipe the centrifuge lid with cleaning solution.
- ➔ Clean the centrifuge chuck and well (avoiding the optical sensor lens) using the disinfectant solution and a lint-free cloth until all traces of blood components are removed.
- ➔ Lubricate the O-ring gasket with a small amount of the silicon lubricant.

Haemonetics recommends that the operator wear protective gloves to avoid direct contact with the cleaning solution and/or any spilled blood which may be present.

In the case of a larger spill, fluid and/or blood may be evacuated into the biohazard waste bag. The operator should complete the following additional steps and contact the local Haemonetics representative for further instructions before using the device:

- ➔ Irrigate the centrifuge drain holes with cleaning solution, until the drain tube is rinsed clear of the spilled material.
- ➔ Remove the bag and replace it with a new bag.
- ➔ Dispose of the used waste bag according to local established policies concerning the disposal of biohazard waste.



Note: A 50 ml syringe of attached to a 20 cm section of disposable tubing placed in the drain holes can be used for irrigation. The biohazard waste bag should be monitored to avoid overfilling.



Warning: An authorized Haemonetics technician should perform a leakage current control after any major fluid spill involving the PCS2 device. Leakage current represents a primary indication of electrical shock hazard and should be checked according to guidelines as described in local standard operating procedures.

Pumps

The pump rotors should be removed from the well with the hexagonal head wrench. Debris should be removed from the rotors and the pump wells on a routine basis, as well as after any spills to contribute to efficient PCS2 operation.

For routine cleaning, the operator should:

- Remove the pump rotor from the housing, using the hexagonal head wrench to remove the pump screw.
- Wipe the rotor and remove all debris from the rollers, using warm water and lint-free cloth or gauze.
- Dry with lint free cloth (or compressed air, if available).
- Clean and dry the pump well using the same method.
- Ensure that all of the rollers spin freely and replace the pump rotor in the well, aligning the cross pin in the rotor with the pump shaft.
- Replace and tighten the hexagonal head screw.

In the case of a fluid spill, the same cleaning method should be followed; however, disinfectant cleaning solution should be used, followed by a clear water rinse. The pump head should not be immersed in water.

Filter screens

The PCS2 device is equipped with filter screens on the bottom of the cabinet, which eliminate dust from incoming cool air. The filters should be cleaned routinely, especially if dust becomes visible on the screens.

To clean the filters, the operator should:

- Remove the retainer plates using a Phillips-head screwdriver.
- Remove the filter screens from the panel.
- Rinse the screens under running water - DO NOT use any cleaning agents.
- Gently squeeze the screens to remove excess water.
- Place the screens on a clean, dry cloth and allow to dry completely.
- Reinsert the screens into the panel, ensuring that all openings are completely covered by the filter.
- Replace the retainer plates and tighten the screws.



Warning: To avoid electrical shock, the filter screen should be completely dry before it is reinstalled on the PCS2 cabinet.

Barcode reader

The barcode reader window should be wiped using a lint-free cloth or gauze and water, then dried. It should be cleaned whenever there is an accumulation of dust or spilled fluid. For optimal cleaning, the operator should remove the protector around the window and replace it once the window has been cleaned.

CUSTOMER SERVICE

Clinical training

Haemonetics employs a staff of Clinical Specialists to provide training for apheresis personnel concerning the use of the PCS2 equipment. The local Haemonetics representative will schedule staff training upon delivery of PCS2 equipment and should be contacted to organize further instruction when needed.

Field service

Haemonetics maintains a worldwide network of company-trained service representatives responsible for responding to technical needs concerning equipment. These technical specialists are available to diagnose and repair any malfunctions, as well as provide routine annual or semi-annual maintenance of the apheresis equipment, including leakage current tests. If service beyond the routine maintenance and cleaning described in this manual is required, the local Haemonetics representative should be contacted to provide specific instruction.

Returned Goods Authorization system

Haemonetics seeks to provide the apheresis customer with equipment and material which respects the highest established standards of quality in design and manufacturing. If for any reason merchandise must be returned to the company, the customer should refer to the Haemonetics Returned Goods Authorization (RGA) system procedure to ensure proper handling and subsequent analysis of the material.

The customer should contact the local Haemonetics Representative (or the Haemonetics Customer Service Department) and provide the following information:

- Product list number, lot number and manufacture date.
- Number of articles to be returned.
- Description of defect.
- Number of parcels being shipped.

The Haemonetics Representative may ask for additional details, depending on the nature of the problem. The customer should be prepared to provide a thorough description of the problem encountered, as well as the product information listed above.

If a contaminated disposable set must be returned by courier services, the Haemonetics Representative may provide specific instructions concerning preparation for shipping blood-contaminated products. In addition to the Haemonetics guidelines, the consumer should strictly follow the local standard operating procedure related to the shipment of blood-contaminated materials and thus minimize any potential health hazards involved.

In some cases, it may be necessary to dispose of the contaminated goods after reporting the problem to the Haemonetics representative. This should be done according to the locally established guidelines pertaining to the disposal of biologically contaminated material.



Warning: Haemonetics products must be properly cleaned and packed prior to their return. It remains an important responsibility of the customer to reduce this serious potential health hazard, by being aware of the risks involved in the shipping, handling and testing of this material.

HAEMONETICS® Cleaning and Maintenance Record

Device serial number: Year of operation:.....

Action																					
Clean cabinet, control panel and valves.																					
Clean air detectors.																					
Clean pump rotors.																					
Clean centrifuge cover and well.																					
Clean optical bowl sensor.																					
Clean optical line sensor.																					
Clean air filters.																					
Inspect O-ring; apply silicon lubricant.																					
Clean DPM (and SPM).																					
Verify biohazard waste bag.																					
Maintenance performed by: (date and initials)																					
Reviewed by: (date and initials)																					

Annual preventive maintenance should be scheduled by a supervisor when appropriate and performed by a Haemonetics service representative or a qualified engineer.

Annual maintenance performed by:..... Date performed:

Reviewed by: Date/initials:.....

Chapter 4

Ensuring Safety and Quality for a PCS2 Procedure

HANDLING THE PCS2 EQUIPMENT	4-2
Storing the PCS2 device and material	4-2
Inspecting the material	4-2
PREVENTING PROBLEMS DURING A PCS2 PROCEDURE	4-3
Understanding the risk of hemolysis	4-3
Avoiding the consequences of flow restriction	4-3
Avoiding bowl misalignment	4-4
Avoiding overheating due to mechanical situations	4-4
Controlling for Red Cell Overrun	4-4
WARNINGS FOR THE OPERATOR	4-5
Electrical shock hazards	4-5
Leakage current control	4-5
Mechanical hazards/rotating parts	4-5
Power outlet connection	4-5
Communicable disease precautions	4-6

HANDLING THE PCS2 EQUIPMENT

Safe and successful PCS2 operation will depend in part on the proper routine handling of the PCS2 equipment. The operator should be aware of the problems which could result if the device or disposable material is stored, installed or used incorrectly.

Storing the PCS2 device and material

The PCS2 device must not be operated or stored in an area where flammable gases or vapors are present. The PCS2 disposable set material should be kept in a dry, well-ventilated area and isolated from any chemical vapors. The operator should handle the disposable set elements with clean, dry hands or gloves.

The ranges for storing the material should be within 8% to 80% rh and -20° C to +50° C. The recommended conditions for the working environment are an ambient temperature between 18° C to 27° C.

Inspecting the material

Prior to installation, the operator should complete a visual inspection of the disposable set elements and control for twisted or flattened sections.

After installing the disposable set, the operator should verify the correct placement of the individual elements, prior to initiating a collection procedure. It is important that the tubing remain free of any twists or occlusions which could cause a flow obstruction.

PREVENTING PROBLEMS DURING A PCS2 PROCEDURE

Understanding the risk of hemolysis

Hemolysis involves the destruction of red blood cell membranes, with the release of free hemoglobin into the plasma portion of the blood. Free hemoglobin does not have the capacity to transport oxygen and can produce serious problems. The remnants of the red cell can stimulate clot formation and damage the vascular nature of the lungs and the kidneys. This could lead to respiratory complications and/or renal failure.

Hemolysis of red cells can occur during an apheresis procedure in the rare event of a mechanically induced situation, such as overheating or excessive pressure.



Warning: Forcing a pump to work against a severe flow restriction can lead to hemolysis, and thus, consequently high levels of free hemoglobin in the plasma. It is important that the operator remain attentive to this fact in the case of any “high return pressure” alarms during PCS2 operation.

If there is any suspicion that hemolysis has occurred, the operator should not return the contents of the bowl to the donor. The local Haemonetics representative should be informed of the problem to provide the operator with further instruction.

Avoiding the consequences of flow restriction

During the DRAW mode, a flow restriction in the effluent tubing can create pressure on the outlet port of the disposable bowl. This unrelieved pressure can deform the rotary seal of the disposable bowl. If the functional characteristics of the rotary seal are altered, the increased friction and excessive heat can make the contents of the bowl unsuitable for return to the donor.

During the RETURN mode, a flow restriction in the effluent tubing can cause the pressure in the centrifuge bowl to drop severely. This sudden drop in pressure could potentially produce hemolysis.

To eliminate these potential problems, the operator should:

- Ensure against inadvertent clamping of the effluent tubing.
- Observe the following changes in flow rate which are possible indications of a flow restriction:
 - Decreased donor flow rate during a Draw cycle.
 - Abnormally increased time required to return the disposable bowl contents bowl to the donor.

Avoiding bowl misalignment



An improperly installed disposable bowl can become misaligned as it spins. This can create excessive friction, and consequently overheat the bowl contents. The operator should verify the alignment of the bowl at the time of installation.

Warning: The operator must not use any bowl which cannot be properly seated in the centrifuge chuck. Overheating can occur, subsequently lead to hemolysis and make any blood being processed unsafe for reinfusion. During operation the operator should interrupt the PCS2 procedure if an abnormality or noise appears, related to the spinning bowl.

Avoiding overheating due to mechanical situations



Overheating could also result from a mechanical or maintenance-related problem, such as a defective bearing or seal within the centrifuge well. In this case, the operator should contact the local Haemonetics representative and discontinue use of the PCS2 device until it is serviced.

Warning: If any component of the PCS2 device has become overheated during a procedure, and thus overheats the blood being processed, the blood components cannot be considered safe for re-infusion.

Controlling for Red Cell Overrun

Red Cell Overrun is the term used to describe the presence of erythrocytes in the effluent tubing and/or product collection container during a plasmapheresis procedure. It is important that the operator observe the appearance of the plasma as it is collected. A pink or reddish hue could indicate a possible red cell spillage or a PCS2 malfunction resulting in severe hemolysis. The situation should be investigated immediately.

If the cause of the reddish hue cannot be determined, the procedure should be discontinued immediately and the blood components in the bowl must not be returned to the donor.

WARNINGS FOR THE OPERATOR



Electrical shock hazards

The operator should always use the PCS2 device with clean, dry hands, or gloves. The internal parts of the PCS2 device contain various electrical components. Contact with any of these components, when the device is connected to an external powered source, could result in an electrical shock to the operator and/or donor.

The operator should never remove any of the PCS2 cabinet panels. Maintenance requiring access to the inner cabinet remains the responsibility of a Haemonetics-trained technician.

Leakage current control

Each PCS2 device receives a careful inspection for leakage current prior to leaving the factory. Haemonetics recommends that a control be performed for leakage current by an authorized representative as part of the annual preventative maintenance. It remains the responsibility of the apheresis center to ensure that this control is performed.

In the event of any major spill in which fluid may enter the cabinet, or an important voltage surge, the operator is responsible to ensure that a leakage current test is performed before re-using the device. The control is necessary to avoid the risk of electrical shock and should be conducted by an authorized Haemonetics representative.

Mechanical hazards/rotating parts

As with any equipment containing rapidly rotating parts, the potential for severe injury exists if contact is made, or if clothing becomes entangled with the moving parts. The PCS2 device contains a safety feature, designed to prevent the centrifuge from spinning if the system has not been properly secured. However, the operator should respect the usual precautions taken when working with equipment containing rotating mechanical parts.

Power outlet connection

To comply with the IEC 60601-1-2 Standard for Medical Electrical Equipment, general requirements for safety, it is not permitted to power the PCS2 device using a multiple portable socket outlet, or an extension cord not supplied by Haemonetics.

Communicable disease precautions

Despite testing and screening to detect communicable diseases such as hepatitis, syphilis or HIV, the risk remains that the blood being processed may be infected. The operator must take the appropriate precautions when handling blood products and disposing of blood-contaminated material, to ensure personal safety as well as the safety of others who may come in contact with the material.

Proper handling of blood contaminated material

If a leak or blood spill should occur, it should be cleaned immediately. The operator should follow the local standard operating procedure outlining the steps to follow and product(s) to be used for the disinfection of blood-contaminated material.

If any blood-contaminated material must be returned to Haemonetics for further inspection, the operator should consult the RGA Procedure as described in the chapter "Maintaining the PCS2 Equipment".

Proper disposal of biologically contaminated materials

Any PCS2 disposable material used during an apheresis process is considered as biologically contaminated. It must be disposed of according to local standard operating procedure for the removal of such material and should not be mixed with non-biologically contaminated waste.